


Enbio PRO


Process times

All process times are drying included.
Start-to-end, from instrument insertion
to having them ready.

Unwrapped FAST / 134°C	Wrapped Class B / 134°C	Wrapped Class B / 121°C	Wrapped PRION Class B / 134°C
9 min	18 min	30 min	43 min



Load capacity per process



Usable load space



Data archiving

All process data is automatically
written to a USB disk.



Process reports

Dedicated software for process
data reading and printing.



Water supply

Enbio uses distilled water from
an external tank.



Power supply

The Enbio Pro requires a
15amp power outlet.



Electrical

Power supply	230V/50Hz
Installed power	3250 W
Max. power consumption	15 A

Pressure and temperature

Working pressure	2.1 Bar
Max. pressure	2.3 Bar
Max. process temperature	137°C

Miscellaneous

Ingress protection	IP20
Data archiving	USB stick
Water conductivity*	< 15 µS/cm

* demineralised water only



Class B
Vacuum
Pouches+
drying



Medical Certified
Medical device class II b
EN-13060 / 98/42/EEC Directive
FDA Approved



Distributed by:

ICONA

DESIGNED WITH YOU IN MIND

Trade Central Business Park
Unit 12/37 Keilor Park Drive
Keilor Park Vic 3042

Ph: 1300 442 662
sales@icona.net.au
www.icona.net.au



Enbio Group AG
Eichengasse 3, 4702 Oensingen
Switzerland

enbio.com